



## **SUTURE PAD 3-SLIT and 5-SLIT**

MD-E-C1000E & MD-E-C1001E

MD-E-C1000S & MD-E-C1001S



User Manual

### **Intended Use**


1. Use our suture pad that simulates animal skin for practicing suturing techniques.
  - 10 uses per pad when used with monofilament

## **1. Parts Included**

1. Pad and case

## **2. Tools Recommended**

1. Monofilament suture
2. Hemostats



### 3. Assembly Instructions

**No Assembly Required.**



## **Best Practices**

To keep your models in best working order we recommend the following:

1. Store your model in a cool, dark, and low dust location.
2. For models containing or exposed to water, wipe down excess, leave all connectors open and allow components to dry completely before storing.
3. If applicable, apply Silicone Finishing Powder to reduce tack and dust accumulation.
4. If applicable, flush or wash models with water to remove lubricant and debris.
5. Work on a flat surface to prevent the model from slipping.

## **Warranty Information**

This is a summary of the Med Dimensions warranty policy provided for your convenience. For the complete and most up-to-date warranty information, please refer to the full policy available at <https://www.med-dimensions.com/warranty/>.

Med Dimensions warrants its products to be free from defects in design, operation, workmanship, and materials. Warranty duration varies by type of product as described on our warranty page. This warranty does not cover damage resulting from misuse, abuse, neglect, or unauthorized repair. To initiate a warranty claim, contact Med Dimensions Customer Support at [support@med-dimensions.com](mailto:support@med-dimensions.com) or [585-510-0005](tel:585-510-0005) x1. Proof of purchase is required for warranty service. This warranty supersedes all other warranties, whether expressed or implied, including warranties of merchantability or fitness for a particular purpose. Med Dimensions is not liable for consequential or incidental damages arising from product use or inability to use.